


**AP-238**

Sample Type: **Biological Fluid**

Latest Revision: **November 11, 2020**

Date of Report: **November 11, 2020**



## 1. GENERAL INFORMATION

<b>IUPAC Name:</b>	1-[4-[(E)-cinnamyl]-2,6-dimethyl-piperazin-1-yl]propan-1-one
<b>InChI String:</b>	InChI=1S/C18H26N2O/c1-4-18(21)20-15(2)13-19(14-16(20)3)12-8-11-17-9-6-5-7-10-17/h5-11,15-16H,4,12-14H2,1-3H3/b11-8+
<b>CFR:</b>	Not Scheduled (11/2020)
<b>CAS#</b>	Not Available
<b>Synonyms:</b>	None Available
<b>Source:</b>	NMS Labs – Toxicology Department

**Important Note:** All identifications were made based on evaluation of analytical data (LC-QTOF-MS) in comparison to analysis of acquired reference material.

**Prepared By:** Alex J. Krotulski, PhD; Melissa F. Fogarty, MSFS, D-ABFT-FT; Donna M. Papsun, MS, D-ABFT-FT; and Barry K. Logan, PhD, F-ABFT

## 2. CHEMICAL AND PHYSICAL DATA

### 2.1 CHEMICAL DATA

Form	Chemical Formula	Molecular Weight	Molecular Ion [M <sup>+</sup> ]	Exact Mass [M+H] <sup>+</sup>
Base	C <sub>18</sub> H <sub>26</sub> N <sub>2</sub> O	286.4	286	287.2118

### 3. SAMPLE HISTORY

AP-238 has been identified in at least one toxicology case since November 2020. The geographical and demographical breakdown are below:

<b>Case Type:</b>	Postmortem (n=1)
<b>Geographical Location:</b>	Indiana (n=1)
<b>Biological Sample:</b>	Femoral Blood (n=1)
<b>Date of First Collection:</b>	August 2020
<b>Additional Findings:</b>	Flualprazolam (n=1), 8-Aminoclonazolam (n=1)

### 4. BRIEF DESCRIPTION

AP-238 is classified as a synthetic opioid and member of the cinnamylpiperazine sub-class. Cinnamylpiperazines are structurally distinct from fentanyl and its analogues. Synthetic opioids, including cinnamylpiperazines, have been linked to adverse events and death. AP-238 is the fourth analogue in this series to be reported by NPS Discovery. AP-238 is structurally similar to AP-237 (bucinnazine), 2-methyl AP-237, and *para*-methyl AP-237. AP-237 is an opioid used therapeutically, although not prescribed within the United States. AP-238, 2-methyl AP-237, and *para*-methyl AP-237 are structural isomers, sharing the same formula and parent mass; however, their chemical behaviors and mass fragmentation patterns differ, allowing for differentiation during analytical testing. AP-238, AP-237, 2-methyl AP-237, *para*-methyl AP-237 and other analogues in this sub-class are not explicitly scheduled in the United States.

### 5. ADDITIONAL RESOURCES

[https://www.caymanchem.com/product/31128/ap-238-\(hydrochloride\)](https://www.caymanchem.com/product/31128/ap-238-(hydrochloride))

[https://www.policija.si/apps/nfl\\_response\\_web/0\\_Analytical\\_Reports\\_final/AP-238-ID-2185-20\\_report.pdf](https://www.policija.si/apps/nfl_response_web/0_Analytical_Reports_final/AP-238-ID-2185-20_report.pdf)

## 6. QUALITATIVE DATA

### 6.1 GAS CHROMATOGRAPHY MASS SPECTROMETRY (GC-MS)


**Testing Performed At:** The Center for Forensic Science Research and Education at the Fredric Rieders Family Foundation (Willow Grove, PA)

**Sample Preparation:** Standard diluted in methanol

**Instrument:** Agilent 5975 Series GC/MSD System

**Standard:** Reference material for AP-238 (Batch: 0592812-2) was purchased from Cayman Chemical (Ann Arbor, MI, USA).  
(<https://www.caymanchem.com/product/31128/ap-238-hydrochloride>)


#### EI (70 eV) Mass Spectrum: AP-238



## 6.2 LIQUID CHROMATOGRAPHY QUADRUPOLE TIME-OF-FLIGHT MASS SPECTROMETRY (LC-QTOF-MS)

<b>Testing Performed At:</b>	The Center for Forensic Science Research and Education at the Fredric Rieders Family Foundation (Willow Grove, PA)
<b>Sample Preparation:</b>	No additional preparation - direct analysis of sample extract
<b>Instrument:</b>	Sciex TripleTOF® 5600+, Shimadzu Nexera XR UHPLC
<b>Column:</b>	Phenomenex® Kinetex C18 (50 mm x 3.0 mm, 2.6 µm)
<b>Mobile Phase:</b>	A: Ammonium formate (10 mM, pH 3.0) B: Methanol/acetonitrile (50:50) Flow rate: 0.4 mL/min
<b>Gradient:</b>	Initial: 95A:5B; 5A:95B over 13 min; 95A:5B at 15.5 min
<b>Temperatures:</b>	Autosampler: 15 °C Column Oven: 30 °C Source Heater: 600 °C
<b>Injection Parameters:</b>	Injection Volume: 10 µL
<b>QTOF Parameters:</b>	TOF MS Scan Range: 100-510 Da Precursor Isolation: SWATH® acquisition (27 windows) Fragmentation: Collision Energy Spread (35±15 eV) MS/MS Scan Range: 50-510 Da
<b>Retention Time:</b>	5.43 min
<b>Standard Comparison:</b>	Reference material for AP-238 (Batch: 0592812-2) was purchased from Cayman Chemical (Ann Arbor, MI, USA). Analysis of this standard resulted in positive identification of the analyte in the extract as AP-238, based on retention time (5.53 min) and mass spectral data. ( <a href="https://www.caymanchem.com/product/31128/ap-238-(hydrochloride)">https://www.caymanchem.com/product/31128/ap-238-(hydrochloride)</a> )

### Extracted Ion Chromatogram (XIC): AP-238



**TOF MS (Top) and MS/MS (Bottom) Spectra: AP-238**

